



DNA-Free Cotton Tipped Applicator

Standard tip, Wood handle



Description

This absorbent cotton tipped applicator is produced with High Quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch shaft by an aqueous based adhesive.

Tested by an independent lab to the most stringent standards (less than 25 pg of human DNA), certified to be DNA-free, DNase-free and RNase-free

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.188" Dia. (4.78 mm)
- 0.625" Length (15.88 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.812" Length (147.62 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip
- Made in USA

Applications

- Forensics / DNA testing
- Crime scene investigation
- Buccal cell collection
- Specimen collection



Puritan®

Ref #	Packaging
25-806 2WC TT FDNA Sterile	10/50/2 (10 boxes per case, 50 tubes per box, 2 swab per tube)