




# HydraFlock® Flocked Swab

Elongated Tip with Polystyrene Handle - no breakpoint



## Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Patented: US, Canada, Europe, Australia, Japan, China  
Patent pending: Hong Kong

## Dimensions

5.93" Overall Length (150.62 mm)

### Tip

- 0.219" Dia. (5.56 mm)
- 0.770" Length (19.56 mm)

### Handle

- 0.100" Dia. (2.54 mm)
- 5.89" Length (149.61 mm)

## Qualities

- Highly absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- Buccal cell collection / DNA testing
- Forensic evidence collection
- Rapid diagnostics
- Specimen collection

| Ref #                 | Packaging   |
|-----------------------|---|
| 3206-H<br>Non-Sterile | 10/100<br>(10 bags per case, 100 per bag)                     |
| 25-3206-H<br>Sterile  | 10/50/1<br>(10 boxes per case, 50 individually wrapped swabs) |



**HydraFlock®**