




# PurFlock Ultra® Flocked Swab

Elongated Tip with Polystyrene Handle - no breakpoint



## Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Patented: US, Canada, Europe, Australia, Japan, China  
Patent pending: Hong Kong

## Dimensions

5.93" Overall Length (150.62 mm)

### Tip

- 0.219" Dia. (5.56 mm)
- 0.770" Length (19.56 mm)

### Handle

- 0.100" Dia. (2.54 mm)
- 5.89" Length (149.61 mm)

## Qualities

- Highly absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- Buccal cell collection / DNA testing
- Forensic evidence collection
- Rapid diagnostics
- Specimen collection

Ref #	Packaging
3206-U Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3206-U Sterile	10/50/1 (10 boxes per case, 50 individually wrapped swabs, 1 swab p/pkg)

