



PurFlock Ultra® Flocked Swab

Elongated Tip with Polystyrene Handle - 20mm breakpoint



Description

An ideal choice in flocked swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

5.93" Overall Length (150.62 mm)

Tip

- 0.219" Dia. (5.56 mm)
- 0.770" Length (19.56 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.89" Length (149.61 mm)

Qualities

- Highly absorbent tip
- Soft smooth tip
- 20mm breakpoint on handle
- Made in USA

Applications

- Buccal cell collection / DNA testing
- Forensic evidence collection
- Rapid diagnostics
- Specimen collection

Ref #	Packaging
25-3206-U 20MM Sterile	10/50/1 (10 boxes per case, 50 individually wrapped swabs, 1 swab p/pkg)

