




HydraFlock® Flocked Swab

Standard Tip with Polystyrene Handle - 80mm breakpoint



Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

6.04" Overall Length (153.416 mm)

Tip

- 0.200" Dia. (5.08 mm)
- 0.69" Length (17.526 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.066" (1.676 mm)
- 6" Length (152.654 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For collection of bacteria or micro-organisms
- Rapid diagnostic testing
- Specimen collection
- Buccal cell collection
- Forensic evidence collection / DNA testing

Ref #	Packaging
3306-H Non-Sterile	10/100 (10 bags of 100 swabs per bag)
25-3306-H Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)



HydraFlock®