




PurFlock Ultra® Flocked Swab

Standard Tip with Polystyrene Handle - 80mm breakpoint



Description

An ideal choice in flock swabs with excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

6.04" Overall Length (153.416 mm)

Tip

- 0.200" Dia. (5.08 mm)
- 0.69" Length (17.526 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.066" (1.676 mm)
- 6" Length (152.654 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For collection of bacteria or micro-organisms
- Rapid diagnostic testing
- Specimen collection
- Buccal cell collection
- Forensic evidence collection / DNA testing

Ref #	Packaging
3306-U Non-Sterile	10/100 (10 bags of 100 swabs per bag)
25-3306-U Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)

