



# Environmental Sampling Swab

EnviroMax® Round dry foam tip, Polypropylene handle



## Description

Ready-to-use environmental foam sampling swabs, provided in wet and dry format for accuracy and consistency of sample collection.

## Dimensions

**6.02" Overall Length (152.908 mm)**

### Tip

- 0.612" Dia. (15.545 mm)
- 2.220" Length (56.388 mm)

### Handle

- 0.125" Dia. (3.175 mm)
- 5.60" Length (142.240 mm)

## Qualities

- Large tip for maximized collection and elution
- Quick-turn-open leak-proof cap
- Swab securely attached to cap
- Rigid foam paddle tip
- High-quality Puritan foam-tipped swab
- Shelf-life: 12 months
- Made in USA

## Applications

- For environmental surface sampling
- Perfect for spot checks
- Quality control
- Detection of microorganisms

Ref #	Packaging
88010-PF RND Non-Sterile	25/1 (25 individually wrapped tubes)
88010-PFB RND Non-Sterile	10/25/1 (10 boxes, 25 individually wrapped tubes, 1 swab per tube)
25-88010 PF RND Sterile	25/1 (25 individually wrapped tubes)
25-88010 PFB RND Sterile	10/25/1 (10 boxes, 25 individually wrapped tubes, 1 swab per tube)

