




HydraFlock® Flocked Swab

Micro Ultrafine Tip with Polystyrene Handle - 100mm breakpoint



HydraFlock®

Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

6.03" Overall Length (153.162 mm)

Tip

- 0.140" Dia. (3.356 mm)
- 0.312" Length (7.925 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.036" (0.914 mm)
- 5.99" Length (152.146 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For nasopharyngeal and urethral cell collection
- Ideal for pediatrics
- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection

Ref #	Packaging
3318-H Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3318-H Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)