



PurFlock Ultra® Flocked Swab

Micro Ultrafine Tip with Polystyrene Handle - 100mm breakpoint



Description

PurFlock's purity make it ideal for PCR, molecular assays, DNA testing, forensic applications as well as direct antigen testing. Its uniquely designed fiber structure accounts for better cell yield in the specimen where sample size is limited.

Patented: US, Canada, Europe, Australia, Japan, China

Patent pending: Hong Kong

Dimensions

6.03" Overall Length (153.162 mm)

Tip

- 0.140" Dia. (3.356 mm)
- 0.312" Length (7.925 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.036" (.914 mm)
- 5.99" Length (152.146 mm)

Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For nasopharyngeal and urethral cell collection
- Ideal for pediatrics
- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection

Ref #	Packaging
3318-U Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3318-U Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)