



HydraFlock® Flocked Swab

Large Tip with Polystyrene Handle - 80mm breakpoint



Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

5.93" Overall Length (150.622 mm)

Tip

- 0.219" Dia. (5.563 mm)
- 0.68" Length (17.27 mm)

Handle

- 0.100" Dia. (2.54 mm) tapered to 0.09" (2.286 mm)
- 5.89" Length (149.606 mm)

Qualities

- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For collection of bacteria or micro-organisms
- Specimen collection
- Buccal cell collection
- Forensic evidence collection

Ref #	Packaging
3406-H Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3406-H Sterile	10/50/1 (10 boxes per case, 50 individually wrapped swabs, 1 swab p/pkg)



HydraFlock®