




PurFlock Ultra® Flocked Swab

Large Tip with Polystyrene Handle - 80mm breakpoint



Description

PurFlock's purity make it ideal for PCR, molecular assays, DNA testing, forensic applications as well as direct antigen testing. Its uniquely designed fiber structure accounts for better cell yield in the specimen where sample size is limited.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

5.93" Overall Length (150.622 mm)

Tip

- 0.219" Dia. (5.563 mm)
- 0.68" Length (17.27 mm)

Handle

- 0.100" Dia. (2.54 mm) tapered to 0.09" (2.286 mm)
- 5.89" Length (149.606 mm)

Qualities

- Highly absorbent tip
- Soft smooth tip
- Made in USA

Applications

- For collection of bacteria or micro-organisms
- Specimen collection
- Buccal cell collection
- Forensic evidence collection

| Ref # | Packaging |
|-----------------------|---|
| 3406-U Non-Sterile | 10/100 (10 bags per case, 100 per bag) |
| 25-3406-U Sterile | 10/50/1 (10 boxes per case, 50 individually wrapped swabs, 1 swab p/pkg) |

