



# HydraFlock® Flocked Swab

Elongated Tip with Polystyrene Handle - 80 mm breakpoint



## Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Patented: US, Canada, Europe, Australia, Japan, China  
Patent pending: Hong Kong

## Dimensions

**5.93" Overall Length (150.622 mm)**

### Tip

- 0.219" Dia. (5.563 mm)
- 0.770" Length (19.558 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 5.89" Length (149.606 mm)

## Qualities

- Medical grade quality
- Highly absorbent tip
- Soft smooth tip
- Made in USA

## Applications

- For Buccal cell collection
- Specimen collection
- Rapid diagnostic testing
- Forensic evidence collection
- DNA testing

Ref #	Packaging
3506-H Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3506-H Sterile	10/50/1 (10 boxes of 50 individually wrapped swabs)

