




HydraFlock® Flocked Swab

Elongated Tip with Polystyrene Handle - no breakpoint



Description

A new era in flock swab technology. HydraFlock provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.219" Dia. (5.56 mm)
- 0.685" Length (17.399 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.067 (1.702 mm)
- 5.96" Length (151.384 mm)

Qualities

- Highly absorbent tip
- Soft smooth tip
- Made in USA
- Handle offered without a breakpoint. Breakpoint can be customized.

Applications

- Buccal cell / DNA collection
- Forensic evidence collection
- Rapid diagnostic testing
- Specimen collection

Ref #	Packaging
3806-H Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3806-H Sterile	10/50/1 (10 boxes per case, 50 individually wrapped swabs, 1 swab p/pkg)



HydraFlock®