




PurFlock Ultra® Flocked Swab

Elongated Tip with Polystyrene Handle - no breakpoint



Description

PurFlock's purity make it ideal for PCR, molecular assays, DNA testing, forensic applications as well as direct antigen testing. Its uniquely designed fiber structure accounts for better cell yield in the specimen where sample size is limited.

Patented: US, Canada, Europe, Australia, Japan, China
Patent pending: Hong Kong

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.219" Dia. (5.56 mm)
- 0.685" Length (17.399 mm)

Handle

- 0.1" Dia. (2.54 mm) tapered to 0.067 (1.702 mm)
- 5.96" Length (151.384 mm)

Qualities

- Highly absorbent tip
- Soft smooth tip
- Made in USA
- Handle offered without a breakpoint. Breakpoint can be customized.

Applications

- Buccal cell / DNA collection
- Forensic evidence collection
- Rapid diagnostic testing
- Specimen collection

Ref #	Packaging
3806-U Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-3806-U Sterile	10/50/1 (10 boxes per case, 50 individually wrapped swabs, 1 swab p/pkg)



**PurFlock
ULTRA®**