



# Foam Tipped Applicator

Foam-over-cotton tip, polypropylene handle



## Description

This absorbent foam tipped applicator is produced with reticulated 100ppi medical grade foam over high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is thermally bonded to a polypropylene handle ensuring the sponge-like foam tip cleans up unwanted particles without releasing any into your critical environment.

## Dimensions

6" Overall Length (152.4 mm)

### Tip

- 0.375" Dia. (9.525 mm)
- 0.937" Length (23.80 mm)

### Handle

- 0.100" Dia. (2.54 mm)
- 6.00" Length (152.40 mm)

## Qualities

- Lint and residue free
- Foam over cotton provides extra absorbency
- Can be used with IPA and DI water, or dry
- Made in USA

## Applications

- General cleaning swab
- Removing contamination
- Applying adhesives and other solutions

Ref #	Packaging
1806-PCF	20/50 (20 resealable bags per case, 50 swabs per bag)



**Puritan®**