



Foam Tipped Applicator

Large foam-over-cotton tip, wooden handle



Description

This absorbent foam tipped applicator is produced with reticulated 100ppi medical grade foam over high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a rigid white birch shaft.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.375" Dia. (9.52 mm)
- 0.999" Length (25.38 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.99" Length (152.28 mm)

Qualities

- Lint and residue free
- Cotton under foam provides extra absorbancy
- Soft smooth tip
- Made in USA

Applications

- General cleaning swab

Ref #	Packaging
1806-WCFL	20/50 (20 bags per case, 50 swabs per bag)
1806-WCFLB	4/250 Bulk (4 bags per case, 250 swabs per bag)



Puritan®