



Foam Tipped Applicator

Polypropylene handle



Description

This small foam tipped applicator is produced with reticulated 100ppi foam. It is heat staked to a black polypropylene handle.

Dimensions

2.70" Overall Length (68.58 mm)

Tip

- 0.156" Dia. (3.962 mm)
- 0.600" Length (15.24 mm)

Handle

- 0.090" Dia. (2.286 mm)
- 2.625" Length (66.675 mm)

Qualities

- Lint and residue free
- Good water absorption
- Soft, non-abrasive paddle tip
- No contaminating residues
- Can be used with IPA, DI water, or dry
- Made in USA

Applications

- Applying adhesives or solutions
- Removing residues and particulates
- General cleaning

Ref #	Packaging
1273 BK	2/500 (2 poly bags per case, 500 swabs per poly bag)



Puritan®