



Foam Tipped Applicator

Standard foam tip, polypropylene handle



Description

This foam tipped applicator is produced with reticulated 100ppi foam. It is thermally bonded to a rigid glass filled polypropylene shaft.

Dimensions

2.375" Overall Length (60.3 mm)

Tip

- 0.188" Dia. (4.77 mm)
- 0.594" Length (15.08 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 2.372" Length (60.26 mm)

Qualities

- Lint and residue free
- Soft and non-abrasive
- Can be used with IPA and DI water, or dry
- Made in USA

Applications

- General cleaning in critical environments
- Removing contamination and residues
- Keyboard cleaning swab

Ref #	Packaging
1802-PF SOLID	20/50 (20 bags per case, 50 swabs per bag)



Puritan®