



Foam Tipped Applicator

Large round foam tip, polypropylene handle



Description

This foam tipped applicator is made with 100ppi medical grade foam. The round head offers a large surface area for effective collection of small samples. This applicator is ideal for a variety of clinical and DNA applications.

Dimensions

5.125" Overall Length (130.175 mm)

Tip

- 0.906" Dia. (23.012 mm)
- 1.125" Length (28.575 mm)

Handle

- 0.125" Dia. (3.175 mm)
- 5.00" Length (127.0 mm)

Qualities

- Round lollipop shape
- Good water absorption
- Lint and residue free
- Made in USA

Applications

- Removing residues or particulates
- Can be used with IPA and DI water, or dry
- General cleaning swab

Ref #	Packaging
1805-PF RND	20/50 (20 resealable bags per case, 50 swabs per bag)



Puritan®