



HydraFlock® Flocked Swab - Dry Transport System

Micro ultrafine tip, polystyrene handle with breakpoint at 100 mm



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. The patented HydraFlock® flocked swab provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.04" (153.416 mm)

Tip

- 0.14" Dia. (3.556 mm)
- 0.312" Length (7.925 mm)

Handle

- 0.100" Dia. (2.54 mm) tapered to 0.036" (0.914 mm)
- 5.99" Length (147.62 mm)

Qualities

- Medical grade quality
- Highly absorbent capacity
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Applications

- For nasopharyngeal and urethral cell collection
- Rapid diagnostic testing
- Forensic evidence collection
- Specimen collection

Ref #	Packaging
25-3318-H BT Sterile	5/100/1 (5 boxes, 100 tubes per box, 1 swab per tube)



Puritan®