



Foam Tipped Applicator

Small rectangular foam tip, polypropylene handle



Description

This foam tipped applicator is produced with reticulated 100ppi Medical Grade foam. It is thermally bonded to a rigid rectangular polypropylene shaft.

Dimensions

5.188" Overall Length (131.78 mm)

Tip

- 0.531" Dia. (13.49 mm)
- 1.00" Length (25.4 mm)

Handle

- 0.125" Dia. (3.18 mm)
- 0.195" Width (4.95 mm)
- 5.00" Length (127 mm)

Qualities

- Lint and residue free
- Soft and non-abrasive
- Can be used with IPA and DI water, or dry
- Made in USA

Applications

- Buccal cell collection
- Applying solutions or solvents
- General cleaning

Ref #	Packaging
1605-PSF RECT Non-Sterile	20/50 (20 bags per case, 50 swabs per bag)
25-1605 1PSF RECT Sterile	10/50/1 (10 boxes of 50 individually wrapped packages, 1 p/pkg)



Puritan®