



# Foam Tipped Applicator

Polypropylene handle scored thumb stop



- Large foam tip

- Scored thumb stop for easy detaching

## Description

This foam tipped applicator is produced with reticulated 100ppi medical grade foam. It is thermally bonded to a rigid rectangular polypropylene handle. Handle has a prominent thumb stop for ease of use when sampling. Scored for easy detaching.

## Dimensions

5.15" Overall Length (130.8 mm)

### Tip

- 0.438" Dia. (11.125 mm)
- 1" Length (25.4 mm)

### Handle

- 0.118" Dia. (2.997 mm)
- 5.025" Length (127.635 mm)

## Qualities

- Medical grade quality
- Prominent thumb stop, scored at 4 cm from tip
- Large foam tip for maximum collection
- Rigid polypropylene handle with thumb stop
- Made in USA

## Applications

- Environmental sample collection for microbial detection and decontamination
- Validated in swab sampling protocol for the detection of bacillus anthracis

Ref #	Packaging
25-1607 1PF SC Sterile	10/50/1 (10 boxes of 50 individually wrapped packages, 1 p/pkg)



**Puritan®**