



Foam Tipped Applicator

Polypropylene handle



- Round foam tip

Description

This foam tipped applicator is produced with reticulated 100ppi medical grade foam. It is thermally bonded to a rigid rectangular polypropylene handle. The round, foam-tipped head offers a large surface area for effective collection of small samples. Ideal for buccal cell collection and luminol test of minute blood stains.

Dimensions

5.125" Overall Length (130.175 mm)

Tip

- 0.906" Dia. (23.012 mm)
- 1.125" Length (28.575 mm)

Handle

- 0.125" Dia. (3.175 mm)
- 5" Length (127 mm)

Qualities

- Medical grade quality
- Large round tip
- High article collection capacity
- Made in USA

Applications

- Specimen collection for DNA testing
- Buccal swab
- Luminol test of minute blood stains

Ref #	Packaging
25-1805 1PF RND Sterile	10/50/1 (10 boxes of 50 individually wrapped packages, 1 p/pkg)



Puritan®