



Foam Tipped Wound Measuring Device

DM Stick, foam tip indented at 0.5 and 1.0 cm, polystyrene handle



Description

This soft foam tipped wound measuring device (DM Stick) features a foam tip on a flexible polystyrene handle. The handle is printed with a graduated metric scale, indented at 1/2 cm intervals. A new and improved foam tip is slightly longer than the original. Tip diameter remains the same, but will now be secured farther along the semi-rigid plastic shaft, adding more soft material at the tip as well as security to the bond.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.2" Dia. (5.08 mm)
- 0.393" Length (9.98 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Soft foam tip for patient comfort
- Handle marked with mm and cm (15 cm overall length)
- Easy measuring
- Made in USA

Applications

- Accurate wound measurement
- Measures: pressure ulcers, stomach and surgical wounds
- Effective for podiatry use
- Used in forensics and evidence collection
- Emergency rooms, long-term care facilities, etc.

Ref #	Packaging
1506-PFB DM Non-sterile	20/50 (20 bags of 50 swabs per bag)
25-1506 1PF DM Sterile	4/50/1 (4 boxes of 50 individually wrapped packages, 1 p/pkg)



Puritan®