



Knitted Polyester Tipped Applicators

Large Paddle Tip, Polypropylene Handle



Description

This lint-free swab has a large knitted polyester fiber tip that is ultrasonically welded to an injection molded polypropylene handle. Ideal for clean room use and other critical environments.

Dimensions

6.185" Overall Length (157.099 mm)

Tip

- 0.39" Dia. (9.906 mm)
- 1.125" Length (28.575 mm)

Handle

- 0.128" Dia. (3.251 mm)
- 6.14" Length (155.956 mm)

Qualities

- Laundered, double ply knitted polyester tip
- Large, flat paddle tip
- Injection molded handle
- Can be used with IPA and DI water or solvents
- Lint-free and non-abrasive
- Made in USA

Applications

- Cleaning with IPA, DI water or other solvents
- Removing residues and particulates
- Applying adhesives or other fluids

Ref #	Packaging
3600 Non-sterile	20/50 (20 bags of 50 swabs per bag = 1,000/ cs)



Puritan®