



HydraFlock® Flocked Swab - Dry Transport System

Large tip, polystyrene handle with breakpoint at 80 mm



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. The patented HydraFlock® flocked swab provides superior collection and release of specimens. The unique microstructure of Hydra's multi-length fibers are designed to enhance rapid absorption and quick release of biological specimens.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 5.93" (150.622 mm)

Tip

- 0.219" Dia. (5.563 mm)
- 0.64" Length (16.256 mm)

Handle

- 0.100" Dia. (2.54 mm) tapered to 0.09" (2.286 mm)
- 5.89" Length (149.606 mm)

Qualities

- Medical grade quality
- Highly absorbent capacity
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Applications

- Buccal cell collection
- Specimen collection
- Forensic evidence collection
- Rapid diagnostic collection

Ref #	Packaging
25-3406-H BT Sterile	5/100/1 (5 boxes, 100 tubes per box, 1 swab per tube)



Puritan®