




Knitted Polyester Tipped Swab

Round Tip, Polypropylene Handle



Description

This lint-free swab has a round, knitted polyester tip that is ultrasonically welded to a Delrin® acetal handle. Perfect for removing epoxy residue in fiber optics manufacturing.

Dimensions

5" Overall Length (127 mm)

Tip

- 0.781" Dia. (19.837 mm)
- 0.203" Length (5.156 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 4.875" Length (123.825 mm)

Qualities

- Lint and residue free
- ATP controlled
- Can be used with IPA, acetone, and other solvents
- No contaminating adhesives
- Made in USA

Applications

- Critical environment cleaning swab
- Used in surface testing for adenosine triphosphate (ATP)
- Suitable for cleanroom environments

Ref #	Packaging
3650 Non-sterile	10/100 (10 resealable poly bags per case, 100 swabs per bag)



Puritan®