




# Cotton Tipped Applicator

Regular tip, Wood handle



## Description

This absorbent double cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a naturally biodegradable wooden handle with a hot-melt glue adhesive.

## Dimensions

6" Overall Length (150.8 mm)

### Tip

- 0.188" Dia. (4.775 mm)
- 0.625" Length (15.875 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.625 mm)

## Qualities

- Medical grade quality
- Regular cotton tip
- Sturdy wood handle, naturally biodegradable
- Hot-melt adhesive not soluble in water
- Made in USA

## Applications

Critical Environments / Industrial

- General purpose swab
- Used for general cleaning

Ref #	Packaging
860-WC Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



**Puritan®**