




Cotton Tipped Applicator

Micro tip, Wood handle



Description

This cotton tipped applicator is produced with high quality USP grade cotton which is carded and coiled onsite to maintain consistency. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. The micro size cotton tip is ideal for hard-to-reach or small areas.

Dimensions

5.875" Overall Length (149.225 mm)

Tip

- 0.125" Dia. (3.175 mm)
- 0.5" Length (12.7 mm)

Handle

- 0.083" Dia. (2.108 mm)
- 5.812" Length (147.625 mm)

Qualities

- Medical grade quality
- Micro size cotton tip for hard-to-reach areas
- Sturdy wood handle, naturally biodegradable
- Lint-free tip
- Made in USA

Applications

Critical Environments / Industrial / Medical / Veterinary

- General, all-purpose swab
- Used for general cleaning of small areas
- Ideal for critical component cleaning
- Used as a podiatry swab

Ref #	Packaging
869-WCS Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



Puritan®