




Cotton Tipped Applicator

Double Small Tip, Fine Paper handle



Description

This double cotton-tipped applicator is produced with high quality USP grade cotton. It is bonded securely to a paper handle by an aqueous based adhesive. It is lint and residue free, making it an ideal applicator for cleaning sensitive components.

Dimensions

2.972" Overall Length (75.489 mm)

Tip

- 0.126" Dia. (3.2 mm)
- 0.551" Length (13.995 mm)

Handle

- 0.06" Dia. (1.524 mm)
- 2.874" Length (73 mm)

Qualities

- Double small tip
- Fine paper handle
- Highly absorbent, does not swell in solvents
- Lint and residue free
- Low in NVR
- Anti-static packaging

Applications

- Cleaning sensitive components
- Cleaning disk drives, optics, computer parts, etc.
- Precision application of solutions or adhesives
- Veterinary
- Weapons cleaning
- Easy disposal

Ref #	Packaging
892-PC DBL Non-sterile	100/25 (100 sleeves per case of 25 swabs per sleeve = 2,500/ cs)



Puritan®



Cotton Tipped Applicator

Double Small Tip, Fine Paper handle

Ref #: 892-PC DBL

Report of Analysis: NVR/Ion/Silicone Oil

Parameter	Result	Unit	Specification
1. Chloride	1.202	µg/ tip	< 8.5
2. NVR (DI)	0.208	mg/ tip	< or = 0.53
3. NVR (IPA)	0.043	mg/ tip	< or = 0.15
4. NVR (Acetone)	0.107	mg/ tip	< or = 0.55
5. Silicone Oil	-	ng/ tip	< 1.3



Puritan[®]