




# Rayon Tipped OB/GYN Applicator

Paper Handle



## Description

Puritan OB/GYN swabs are made with USP grade rayon fiber that is carded on-site to maintain consistency. The large, oversized, rayon tip is extra-absorbent.

8.188" Overall Length (207.975 mm)

### Tip

- 0.5" Dia. (12.7 mm)
- 1.125" Length (28.575 mm)

### Handle

- 0.156" Dia. (3.962 mm)
- 8" Length (203.2 mm)

## Qualities

- Extra absorbent, large oversized tip
- Long paper handle
- Medical grade quality
- Made in USA

## Applications

- Preparatory to specimen sampling or exams
- Forensics/evidence collection
- For cell collection in veterinary applications

Ref #	Packaging
808 Non-sterile	10/50/1 (10 boxes of 50 swabs each, 1 p/pkg)
808 100 Non-sterile	10/50/2 (10 boxes of 50 swabs each, 2 p/pkg)
808 BULK Non-sterile	500 Bulk (10 bundles per case, 50 per bundle)



**Puritan®**