



Rayon Tipped Proctoscopic Applicator

Polypropylene handle



Puritan®

Description

Puritan Proctoscopic swabs are made with USP grade rayon fiber that is carded on-site to maintain consistency. The large rayon tip is extra-absorbent.

Dimensions

16.062" Overall Length (407.975 mm)

Tip

- 0.5" Dia. (12.7 mm)
- 1.125" Length (28.575 mm)

Handle

- 0.150" Dia. (3.81 mm)
- 15.875" Length (403.225 mm)

Qualities

- Extra absorbent, large oversized tip
- Medical grade quality
- Made in USA

Applications

- Designed for use with proctoscopes
- OB/GYN preparation & examinations
- Gynecological cleaning
- Laboratory cell research
- For cell collection in veterinary applications

Ref #	Packaging
816 Non-sterile	10/500 (10 boxes per case, 50 per box)
816 100 Non-sterile	10/100 (10 boxes per case, 100 per box)
816 BULK Non-sterile	500 Bulk (10 bundles per case, 50 per bundle)