




Rayon Tipped Proctoscopic Applicator

Polypropylene handle



Description

Puritan proctoscopic swabs are made with USP grade rayon fiber that is carded on-site to maintain consistency. Extra binder is used for a firmer rayon tip.

Dimensions

16.062" Overall Length (407.975 mm)

Tip

- 0.5" Dia. (12.7 mm)
- 1.125" Length (28.575 mm)

Handle

- 0.150" Dia. (3.81 mm)
- 15.875" Length (403.225 mm)

Qualities

- Firmer tip and more tightly wound than the standard 816 swab
- Multi-use swab
- Medical grade quality
- Made in USA

Applications

- Preparatory to specimen sampling or exams
- Design for use with proctoscopes
- Gynecological cleaning
- Laboratory cell research
- For cell collection in veterinary applications

Ref #	Packaging
817 Non-sterile	10/50 (10 boxes of 50 swabs each)



Puritan®