



# Rayon Tipped OB/GYN & Proctoscopic Applicator

Polystyrene handle



-USP grade rayon tip

\*Image not actual size

## Description

This specialty swab is made with USP grade rayon fiber that is carded on-site to maintain consistency. The applicator has an oversized double tip that is bonded to a scored polystyrene handle. When this 8" double tipped swab is broken at the score location it results in two 4" swabs.

## Dimensions

8" Overall Length (203.20 mm)

### Tip

- 0.437" Dia. (11.1 mm)
- 1.125" Length (28.575 mm)

### Handle

- 0.15" Dia. (3.81 mm)
- 7.75" Length (196.85 mm)

## Qualities

- Oversized double rayon tip
- Scored polystyrene handle
- Medical grade quality
- Made in USA

## Applications

- Used for pre-saturation
- Packaged with antiseptics

Ref #	Packaging
818 Non-sterile	2,000 Bulk



**Puritan**<sup>®</sup>