




# Cotton Tipped OB/GYN Applicator

Paper handle



-USP grade cotton tip

- Paper handle

\*Image not actual size

## Description

This Puritan OB/GYN swab is made with High Quality USP grade cotton which is carded on-site to maintain consistency. The large, oversized, cotton tip is extra-absorbent.

## Dimensions

8.188" Overall Length (207.975 mm)

### Tip

- 0.5" Dia. (12.7 mm)
- 1.125" Length (28.575 mm)

### Handle

- 0.156" Dia. (3.962 mm)
- 8" Length (203.2 mm)

## Qualities

- Extra absorbent, large oversized tip
- Long paper handle
- Medical grade quality
- Made in USA

## Applications

- Preparatory to specimen sampling or exams
- General cleaning of large surfaces

| Ref #                     | Packaging                                 |
|---------------------------|---|
| 808 COTTON<br>Non-sterile | 10/500<br>(10 boxes per case, 50 per box) |



**Puritan®**