



Polyester Tipped Applicator

Miniature tip, ultrafine polystyrene handle



Description

This specialty applicator is tipped with medical grade polyester that we card ourselves. They're non-absorbent and contain no fatty acids or latent sugars, making them ideal for transport, lab work, and diagnostic kits. The miniature tip is ideal for male specimen collection including STD testing.

Dimensions

6.062" Overall Length (153.975 mm)

Tip

- 0.078" Dia. (1.981 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.098" Dia. (2.489 mm) tapered to 0.052" (1.321 mm)
- 6" Length (152.4 mm)

Qualities

- Medical grade quality
- Synthetic polyester contains no fatty acids or latent sugars
- Thin flexible polystyrene handle allows for minimal patient discomfort
- Made in USA

Applications

- Male specimen collection
- STD testing
- Testing and screening swab

Ref #	Packaging
25-800 1PD 50 Sterile	10/50/1 (10 boxes of 50 individually wrapped packages, 1 swab p/pkg)



Puritan®