




Polyester Tipped Applicator

Aluminum/plastic handle



- Mini-tip

- Aluminum/plastic handle

Description

This specialty applicator is tipped with medical grade polyester that we card ourselves. They're non-absorbent and contain no fatty acids or latent sugars, making them ideal for transport, lab work, and diagnostic kits. The miniature polyester tip is adhered to flexible aluminum which is inserted into a rigid plastic handle.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.078" Dia. (1.981 mm)
- 0.688" Length (17.475 mm)

Handle

- 0.098" Dia. (2.489 mm) tapered to 0.035" (0.889 mm)
- 6" Length (152.4 mm)

Qualities

- Medical grade quality
- Synthetic polyester contains no fatty acids or latent sugars
- Thin flexible aluminum on rigid plastic handle
- Made in USA

Applications

- Specimen collection
- Testing and screening swab

Ref #	Packaging
800-PD ALUM Non-Sterile	10/100 (10 bags per case, 100 per bag)
25-800 1PD ALUM 50 Sterile	10/50/1 (10 boxes of 50 individually wrapped packages, 1p/pkg)



Puritan®