



Polyester Tipped Applicator

Mini-tip, Aluminum handle



25-800 D 50



Description

This specialty applicator is tipped with medical grade polyester that we card ourselves. They're non-absorbent and contain no fatty acids or latent sugars, making them ideal for transport, lab work, and diagnostic kits.

Dimensions

5.5" Overall Length (139.7 mm)

Tip

- 0.078" Dia. (1.981 mm)
- 0.687" Length (17.475 mm)

Handle

- 0.035" Dia. (0.889 mm)
- 5.375" Length (136.525 mm)

Qualities

- Medical grade quality
- Synthetic polyester contains no fatty acids or latent sugars
- Thin flexible aluminum handle follows contours of passages with virtually no injury to mucosal tissues
- Made in USA

Applications

- Forensic evidence collection
- Diagnostic specimen collection for testing and screening
- Flu and strep testing
- Pediatric use
- Eye microbiology
- Veterinary care

Ref #	Packaging
25-800 D 50 Sterile	10/50/1 (10 boxes of 50 individually wrapped packages, 1 p/pkg)



Puritan®