



Polyester Tipped Applicator

Aluminum handle, curled end



Description

This specialty applicator is tipped with medical grade polyester that we card ourselves. They're non-absorbent and contain no fatty acids or latent sugars, making them ideal for transport, lab work, and diagnostic kits.

Dimensions

5.437" Overall Length (139.7 mm)

Tip

- 0.078" Dia. (1.981 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.025" Dia. (0.635 mm)
- 5.375" Length (136.525 mm)

Qualities

- Medical grade quality
- Synthetic polyester tip contains no fatty acids or latent sugars
- Thin aluminum handle has curled end which acts as a handle
- Ideal for diagnostic use
- Made in USA

Applications

- Male specimen collection
- Specimen collection for testing and screening
- Applicator suitable for small cavities
- STD testing
- Pediatric ENT
- Eye microbiology
- Veterinary care
- Avian disease monitoring

Ref #	Packaging
25-801 D 50 Sterile	10/50/1 (10 boxes of 50 individually wrapped packages, 1 swab p/pkg)



Puritan®