



Polyester Tipped Applicator

Polystyrene handle



Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

Dimensions

3" Overall Length (76.2 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 2.937" Length (74.599 mm)

Qualities

- Medical grade quality
- Absorbent tip
- Soft smooth tip
- Colored handle for coded uses
- Regular size tip
- Made in USA

Applications

- DNA Collection
- Paternity testing
- Diagnostic testing



Ref #	Packaging
25-803 2PD BLUE Sterile	10/100/2 (10 boxes of 100 individually wrapped packages, 2p/pkg)
25-803 2PD PINK Sterile	
25-803 2PD YELLOW Sterile	