



PurFlock Ultra® Flocked Swab - Dry Transport System

Elongated tip, polystyrene handle with breakpoint at 30 mm



Description

Puritan's polypropylene dry transport tube and plug-style cap is perfect for specimen collection and transportation to the lab or use in diagnostic kit. PurFlock Ultra® flocked swabs provide excellent absorption and elution of specimens when HydraFlock's high levels of absorbency aren't necessary. PurFlock Ultra's high purity characteristics makes it an ideal swab for various rapid diagnostic tests.

Dimensions

Overall Length

- Tube: 6.062" (153.975 mm)
- Swab: 6.00" (152.4 mm)

Tip

- 0.219" Dia. (5.563 mm)
- 0.685" Length (17.399 mm)

Handle

- 0.100" Dia. (2.54 mm) tapered to 0.067" (1.702 mm)
- 5.96" Length (151.384 mm)

Qualities

- Medical grade quality
- Highly absorbent capacity
- Packaged in a reclosable dry transport tube for ease in transport to the lab
- Made in USA

Applications

- Rapid diagnostic testing
- Forensic evidence collection
- Specimen collection
- For buccal cell collection
- DNA testing

Ref #	Packaging
25-3606-U BT Sterile	5/100/1 (5 boxes, 100 tubes per box, 1 swab per tube)



Puritan®