



Polyester Tipped Applicator

Polystyrene handle



- Regular size tip

- Semi-flexible polystyrene handle

* Not actual size. 25-806 2PD has 2 swabs per package.



Puritan®

Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It is non-absorbent and contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive. Available sterile or non-sterile.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.23 mm)

Qualities

- Medical grade quality
- Polyester tip contains no fatty acids or latent sugars
- Regular size tip
- Semi-flexible polystyrene handle
- Made in USA

Applications

- Diagnostic specimen collection
- Patient care
- Cell collection for DNA testing
- Gynecological screening

Ref #	Packaging
806-PD Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-806 1PD Sterile	10/100/1 (10 boxes per case of 100 individually wrapped packages, 1 per/pkg)
25-806 2PD Sterile	10/100/2 (10 boxes per case of 100 individually wrapped packages, 2 per/pkg)