



Polyester Tipped Applicator

Polystyrene handle



- Regular size tip

- Solid polystyrene handle



Puritan®

Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It is non-absorbent and contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a solid polystyrene handle by an aqueous based adhesive. Inherent score mark on handle.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.100" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Low absorbency tip
- Synthetic polyester contains no fatty acids or latent sugars
- Non-migrating handle (no fluid intake in handle)
- Made in USA

Applications

- Diagnostic specimen collection
- Patient care
- Cell collection for DNA testing
- Gynecological screening

Ref #	Packaging
25-806 1PD SOLID Sterile	10/100/1 (10 boxes of 100 individually wrapped-packages, 1 swab p/pkg)
25-806 2PD SOLID Sterile	10/100/2 (10 boxes of 100 individually wrapped-packages, 2 swabs p/pkg)