




# Polyester Tipped Applicator

Wood handle



- Standard size tip

- Sturdy wood handle



**Puritan®**

## Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It is non-absorbent and contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a sturdy wood handle by an aqueous based adhesive. Available sterile or non-sterile.

## Dimensions

6" Overall Length (152.4 mm)

### Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

### Handle

- 0.1" Dia. (2.54 mm)
- 5.812" Length (147.625 mm)

## Qualities

- Medical grade quality
- Polyester tip contains no fatty acids or latent sugars
- Sturdy wood handle, naturally biodegradable
- Made in USA

## Applications

- Specimen collection
- Testing and screening swab
- Patient care
- Surface collection

Ref #	Packaging
806-WD Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-806 1WD Sterile	10/100/1 (10 boxes per case of 100 individually wrapped packages, 1 p/pkg)
25-806 2WD Sterile	10/100/2 (10 boxes per case of 100 individually wrapped packages, 2 p/pkg)