



Polyester Tipped Applicator

Polystyrene handle



- Regular size tip

- Semi-flexible polystyrene handle

Description

This synthetic tipped applicator is made of medical-grade polyester fiber that is carded in our factory to maintain consistency. It contains no fatty acids or latent sugars, making it ideal for specimen transport, lab work, and diagnostic kits. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

Dimensions

3" Overall Length (76.2 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 2.937" Length (74.599 mm)

Qualities

- Medical grade quality
- Polyester tip contains no fatty acids or latent sugars
- Regular size tip
- Made in USA

Applications

- Cell and surface collection
- Patient care

Ref #	Packaging
803-PD Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



Puritan®