



Polyester Tipped Applicator

Polystyrene handle



- Regular size tip

- Semi-flexible polystyrene handle

* Not actual size

Description

This synthetic tipped applicator is made of ATP controlled polyester fiber that is carded in our factory to maintain consistency. It is bonded securely to a semi-flexible polystyrene handle by an aqueous based adhesive.

Dimensions

5" Overall Length (130.175 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.031" Length (127.787 mm)

Qualities

- ATP controlled
- Polyester tip contains no fatty acids or latent sugars
- Regular size tip
- Semi-flexible polystyrene handle
- Made in USA

Applications

- Surface testing for adenosine triphosphate
- Patient care

Ref #	Packaging
805-PD Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)



Puritan[®]