



# Rayon Tipped Applicator

Aluminum handle

800 R

25-800 R 50



- Flexible aluminum handle

## Description

This tipped applicator is made of the finest USP grade rayon fiber that is carded in our factory to maintain consistency. Ideal for diagnostic specimen collection, wound treatment, and application of medicines. Flexible aluminum handle for minimal patient discomfort. Available sterile or non-sterile.

## Dimensions

	Overall Length	Tip (length x dia.)	Handle (length x dia.)
800 R	6.547" 165.1 mm	0.687 x 0.078" 17.45 x 1.981 mm	6.5 x 0.035" 165.1 x 1.981 mm
25-800 R 50	5.5 " 139.7 mm		5.375 x 0.035" 136.525 x 0.889 mm

## Qualities

- Medical grade quality
- Flexible aluminum handle
- Made in USA

## Applications

- Specimen collection
- General patient care
- Wound treatment
- Applying medicines

Ref #	Packaging
800 R Non-sterile	10/100 (10 bags of 100 swabs per bag)
25-800 R 50 Sterile	10/50/1 (10 boxes of 50 individually wrapped packages, 1 p/pkg)



**Puritan®**