




Rayon Tipped Applicator

Aluminum handle



Description

This specialty tipped applicator is made of the finest USP grade of rayon fiber that is carded in our factory to maintain consistency. Ideal for diagnostic specimen collection, wound treatment, and application of medicines. Flexible aluminum handle to minimize patient discomfort.

Dimensions

5.437" Overall Length (138.1 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.025" Dia. (0.635 mm)
- 5.375" Length (136.525 mm)

Qualities

- Medical grade quality
- Flexible aluminum handle
- Made in USA

Applications

- Specimen collection
- General patient care
- Wound treatment
- Applying medicines

Ref #	Packaging
25-801 R 50 Non-sterile	10/50/1 (10 boxes per case, 50 individually wrapped packages, 1 p/pkg)



Puritan®