



Rayon Tipped Applicator

Polystyrene handle



*Not actual size



Puritan®

Description

This specialty swab is made with USP grade rayon fiber that is carded on-site to maintain consistency. It is bonded securely to a rigid white birch handle by an aqueous based adhesive. Available sterile or non-sterile.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.203" Dia. (5.156mm)
- 0.687" Length (17.45 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Semi-flexible polystyrene handle
- Made in USA

Applications

- Specimen collection swab for testing and screening
- General purpose collection swab
- Cell collection for DNA testing

Ref #	Packaging
806-PR Non-sterile	10/10/100 (10 boxes per case, 10 bags per box, 100 swabs per bag)
25-806 1PR Sterile	10/100/1 (10 boxes per case of 100 individually wrapped packages, 1 per/pkg)
25-806 2PR Sterile	10/100/2 (10 boxes per case of 100 individually wrapped packages, 2 per/pkg)