



Rayon Tipped Applicator

Regular tip, Wood handle



- Regular size tip

*Not actual size.

Description

This specialty swab is made with USP grade rayon fiber that is carded on-site to maintain consistency. This swab is ideal for diagnostic specimen collection, wound treatment and application of medicines. It is bonded securely to a sturdy wood handle by an aqueous based adhesive.

Dimensions

6" Overall Length (152.4 mm)

Tip

- 0.203" Dia. (5.156 mm)
- 0.687" Length (17.45 mm)

Handle

- 0.1" Dia. (2.54 mm)
- 5.875" Length (149.225 mm)

Qualities

- Medical grade quality
- Sturdy wood handle, naturally biodegradable
- Made in USA

Applications

- General purpose applicator
- Specimen collection
- Application of topical dressings and first aid
- Industrial general cleaning
- Cleaning sensors and components
- Removing residues

Ref #	Packaging
25-806 1WR Sterile	10/100/1 (10 boxes per case of 100 individually wrapped packages, 1 p/pkg)



Puritan®